

1. Part Number: 5065C-1001B

2. Photo:



Pendulum Beaker fitted to a Telescopic Rod.
Note: Telescopic Rod is sold separately.



Pendulum Beaker 1000ml
(Blue rotating elbow is included in the price)

3. Product Information:

Description	Pendulum Beaker 1000ml – 316Lss
Materials:	Beaker: 316L stainless steel Rotating Elbow PBT - Blue
Surface Finish:	Better than 0.8 microns Ra
Nominal Volume:	1000ml
Internal Beaker Diameter:	100mm
Overall Beaker Diameter:	110mm
Nominal Beaker Height:	145mm
Nominal Weight:	750g
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: Suitable for use with 5065A and 5065S telescopic rods.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 11-Dec-24

