

1. Part Number: 5065C-1002

2. Photo:



Pendulum Beaker with Telescopic Rod.
Note: Telescopic Rod is sold separately.



Pendulum Beaker
(Beaker + Rotating Elbow)

3. Product Information:

Description	Pendulum Beaker 1000ml	
Materials of Construction:	Beaker	Polypropylene - Translucent
	T-Bar	Polypropylene - Blue
	Rotating Elbow	PBT – Blue
Nominal Volume:	1000ml	
Maximum Beaker Diameter:	115mm	
Overall Beaker Width:	145mm	
Nominal Beaker Height:	148mm	
Nominal Weight:	200g	
Recommended Storage Conditions:	Dry and ambient temperature	

Note: The Rotating Elbow will fit into the shank on the side of the Beaker.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 24-Jun-24