


1. Part Number: 5065F-1000

2. Photograph:



3. Product Information:

Description:	Chemical Dipper 1000ml
Material of Construction:	Polypropylene
Nominal Volume:	1000ml
Surface Finish:	N/A
Diameter of Beaker:	117mm
Length of Handle:	980mm
Overall Length:	1190mm
Nominal Weight:	360g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-Dec-24