


1. Part Number: 5065F-250

2. Photograph:



3. Product Information:

Description:	Chemical Dipper 250ml
Material of Construction:	Polypropylene
Nominal Volume:	250ml
Surface Finish:	N/A
Diameter of Beaker:	79mm
Length of Handle:	980mm
Overall Length:	1130mm
Nominal Weight:	305g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-Dec-24