

**1. Part Number:** 5065F-500

**2. Photograph:**



**3. Product Information:**

Description:	Chemical Dipper 500ml
Material of Construction:	Polypropylene
Nominal Volume:	1000ml
Surface Finish:	N/A
Diameter of Beaker:	90mm
Length of Handle:	980mm
Overall Length:	1170mm
Nominal Weight:	325g
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-Dec-24