

1. Part Number: 5100A-1000

2. Photo:



3. Product Information:

Description	PTFE Dipper 1000ml
Product Contact Parts:	PTFE (Polytetrafluoroethylene) <ul style="list-style-type: none"> • Pure, virgin PTFE • USP Class VI compliant • FDA conforming
Nominal Capacity:	1000ml
Diameter of Cup:	105mm
Handle Length:	600mm
Overall Length:	775mm
Working Temperature Range:	-200 to +280°C
Nominal Weight:	845g
Recommended Storage Conditions:	Dry and ambient temperature

Notes: Increase the effective reach of this sampler by using Extension Rods

- PTFE Extension Rod 600mm Part No. 5100E-600
- PTFE Extension Rod 1000mm Part No. 5100E-1000

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 1-Oct-19