


1. Part Number: 5100A-500**2. Photo:****3. Product Information:**

Description	PTFE Dipper 500ml
Product Contact Parts:	PTFE (Polytetrafluoroethylene) <ul style="list-style-type: none">• Pure, virgin PTFE• USP Class VI compliant• FDA conforming
Nominal Capacity:	500ml
Diameter of Cup:	80mm
Handle Length:	600mm
Overall Length:	740mm
Working Temperature Range:	-200 to +280°C
Nominal Weight:	655g
Recommended Storage Conditions:	Dry and ambient temperature

Notes:

- Increase the effective reach of this sampler by using Extension Rods
- PTFE Extension Rod 600mm Part No. 5100E-600
- PTFE Extension Rod 1000mm Part No. 5100E-1000

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 1-Oct-19