

**1. Part Numbers:** 5132S-300

**2. Photograph:**



**3. Product Information:**

Description	Immersion Cylinder STAINLESS – 300ml
Material of Construction:	304 Stainless Steel
Nominal Volume:	300ml
Body Length	129mm
Body Diameter:	70mm
Nominal Weight:	1.3kg
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 27-Oct-20