


1. Part Number: 7008A-01

2. Photo:



3. Product Information:

Description	MiniCorer 150mm
Materials:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	190mm
Blade Length:	150mm
Max. Blade Width:	19mm
Nominal Weight:	115g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 31-Mar-20