

1. Part Numbers: 8011K-02

2. Photograph:



3. Product Information:

Description	TeleSwab Kit 2
Material of Construction:	Pole = Aluminium Swab Template = 304 stainless steel & PTFE
Minimum Length:	830 mm
Maximum Length:	3200 mm
Number of Sections:	5
Weight:	850 g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-Jan-25