


	<h1>Product Specification</h1>	Specification Number	8011K-02
		Revision	01

Part No. 8011K-02

Description TeleSwab Kit 2

Material of Construction: Pole = Aluminium
Swab Template = 304 stainless steel & PTFE




Minimum Length: 830 mm

Maximum Length: 3200 mm

Number of Sections: 5

Weight: 850 g

Recommended Storage Conditions: Dry and ambient temperature

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 19-Feb-19
Page 1 of 1	Sampling Systems Ltd. 4-5 Forge Mills Park, Station Road, Coleshill, Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com	