

1. Part Number: 8022P-100

2. Photograph:



3. Product Information:

Description:

Liquid Funnel 100mm

Material of Construction:

Polypropylene (PP) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction:

Single piece, injection moulded

Max. Diameter:

117mm

ID at Top of Funnel:

100mm

Overall Height:

103mm

OD at Bottom of Spout:

11mm

ID at Bottom of Spout:

8mm

OD at Top of Spout:

17mm

Nominal Weight:


16g

Recommended Storage Conditions:

Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-Jun-24