

**1. Part Number:** 8022P-350

**2. Photograph:**



### 3. Product Information:

Description:

Liquid Funnel - 350mm

Material of Construction:

Polypropylene (PP) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction:

Single piece, injection moulded

Maximum Diameter:

343mm

Overall Height:

264mm

Outside Diameter of Tip of Spout:

52mm

Inside Diameter of Tip of Spout:

48mm

Length of Spout:

43mm

Nominal Weight:

348g

Recommended Storage Conditions:

Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-Jun-24