

**1. Part Number:** 8022P-50

**2. Photograph:**



**3. Product Information:**

**Description:**

Liquid Funnel 50mm

**Material of Construction:**

Polypropylene (PP) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

**Method of Construction:**

Single piece, injection moulded

**Max. Diameter:**

68mm

**ID at Top of Funnel:**

50mm

**Overall Height:**

56mm

**OD at Bottom of Spout:**

8mm

**ID at Bottom of Spout:**

5mm

**OD at Top of Spout:**

12mm

**Nominal Weight:**

19g

**Recommended Storage Conditions:**

Dry and ambient temperature

**4. Approval:**

Signed:



**Name:** Kuldeep Flora

**Position:** Quality Manager

**Date:** 21-Jun-24