

1. Part Number: 8023H-200

2. Photograph:



3. Product Information:

Description:

Powder Funnel - 200mm

Material of Construction:

High Density Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction:

Single piece, injection moulded

Maximum Diameter:

200mm

Overall Height:

165mm

OD at Bottom of Spout:

43mm

ID at Bottom of Spout:

40mm

Nominal Weight:

86g

Recommended Storage Conditions:

Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 4-Mar-20