

**1. Part Number:** 8023P-100

**2. Photograph:**



**3. Product Information:**

Description:	Powder Funnel - 100mm
Material of Construction:	Polypropylene (PP) - virgin <ul style="list-style-type: none"><li>• Conforms to FDA 21 CFR 177.1520</li><li>• Conforms to EU Regulations 10/2011</li><li>• Conforms to EC Regulations 1935/2004</li><li>• BSE/TSE free</li></ul>
Method of Construction:	Single piece, injection moulded
Max. Top Diameter	109mm
Overall Spout Diameter	20mm
Internal Spout Diameter	17mm
Spout Length	32mm
Overall Height	106mm
Nominal Weight:	29g
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 28-Sep-20