

	<h1 style="margin: 0;">Sampler Specification</h1>	Specification Number	8023P-150
		Revision	01

Part No. 8023P-150

Description Powder Funnel – 150mm Top Diameter

Material of Construction: Polypropylene (PP)



Max. Top Diameter 150mm

Overall Spout Diameter 28mm

Internal Spout Diameter 24mm

Spout Length 30mm


Overall Height 135mm

Individually Bagged? No

Pre-sterilised? No

Weight 59g

Number of Samplers per Box: Sold individually

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 24-Mar-17
----------------------------	--	--