

1. Part Number: 8023P-350

2. Photograph:



3. Product Information:

Description:	Powder Funnel - 350mm
Material of Construction:	Polypropylene (PP) - virgin <ul style="list-style-type: none"> • Conforms to FDA 21 CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Maximum Diameter:	343mm
Overall Height:	264mm
Outside Diameter of Tip of Spout:	52mm
Inside Diameter of Tip of Spout:	48mm
Length of Spout:	43mm
Nominal Weight:	348g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Jonathan Fenner

Position: Director

Date: 3-Nov-22