

**1. Part Number:** 8045A-250S

**2. Photograph:**



**3. Product Information:**

Description:	SteriWare PharmaScoop® 250ml - Sterile
Use:	Single use, disposable scoop
Material of Construction:	High Impact Polystyrene (HIPS) – rigid, white <ul style="list-style-type: none"> <li>• FDA acceptable CFR 177.1640</li> <li>• EU Regulations 10/2011 compliant</li> <li>• EC regulation 1935/2004 compliant</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Overall Length:	209mm
Handle Length:	104mm
Maximum Width:	81mm
Nominal Weight:	45g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Number of Scoops per Box:	85
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

#### 4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

#### 5. Country of Origin:

Made and packed in the UK.

#### 6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

#### 7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-Jul-24