

	<h1>Product Specification</h1>	Specification Number	8045A-30
		Revision	06


Part No. 8045A-30

Description SteriWare PharmaScoop® – 30ml

Use: Single use, disposable scoop

Material of Construction: High Impact Polystyrene (HIPS) – rigid, white

- FDA acceptable CFR 177.1640
- EU Regulations 10/2011 compliant
- EC regulation 1935/2004 compliant



Method of Construction: Single piece, injection moulded

Overall Length: 153mm

Handle Length: 101mm

Maximum Width: 40mm

Nominal Weight: 10g

Assembly & Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE bag)


Method of Sterilisation: N/A

Number of Samplers per Box: 100

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

SteriWare and PharmaScoop are registered Trademarks belonging to Sampling Systems Ltd.

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 3-Oct-17
----------------------------	----------------------------------------------------------------------------------------------------	---------------------------------------------------------------