

**1. Part Number:** 8045B-2500S

**2. Photograph:**



### 3. Product Information:

|                                 |   |
|---------------------------------|---|
| Description:                    | SteriWare® PharmaScoop – 2500ml - Sterile   |
| Use:                            | Single use, disposable scoop  |
| Material of Construction:       | Polypropylene (PP) – virgin – blue <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul> |
| Method of Construction:         | Single piece, injection moulded   |
| Nominal Volume:                 | 2500ml  |
| Overall Length:                 | 357mm   |
| Handle Length:                  | 147mm   |
| Maximum Width:                  | 166mm   |
| Nominal Weight:                 | 224g  |
| Moulding Environment:           | Class 100,000 Cleanroom   |
| Primary Packing Environment:    | Class 10,000 Cleanroom  |
| Individually Bagged?            | Yes (heat sealed PE/PET bag)  |
| Method of Sterilisation:        | Gamma Irradiation (25 - 45 kGy)   |
| Quantity per Box:               | 15  |
| Recommended Storage Conditions: | Dry and ambient temperature   |
| Shelf Life:                     | 5 years from the date of manufacture  |
| ATEX:                           | This device has not been ATEX tested  |

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-Jul-24