

1. Part Number: 8045B-500S

2. Photograph:



3. Product Information:

Description:	SteriWare® PharmaScoop – 500ml - Sterile
Use:	Single use, disposable scoop
Material of Construction:	Polypropylene (PP) – virgin – blue <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	500ml
Overall Length:	249mm
Handle Length:	118mm
Maximum Width:	97mm
Nominal Weight:	86g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25 to 45kGy)
Quantity per Box:	50
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-Jul-24