

1. Part Number: 8045B-500

2. Photograph:



3. Product Information:

| | |
|---------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description: | SteriWare® PharmaScoop – 500ml |
| Use: | Single use, disposable scoop |
| Material of Construction: | Polypropylene (PP) – virgin – blue <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free |
| Method of Construction: | Single piece, injection moulded |
| Nominal Volume: | 500ml |
| Overall Length: | 249mm |
| Handle Length: | 118mm |
| Maximum Width: | 97mm |
| Nominal Weight: | 86g |
| Moulding Environment: | Class 100,000 Cleanroom |
| Primary Packing Environment: | Class 10,000 Cleanroom |
| Individually Bagged? | Yes (heat sealed PE/PET bag) |
| Method of Sterilisation: | N/A |
| Quantity per Box: | 50 |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Shelf Life: | 5 years from the date of manufacture |
| ATEX: | This device has not been ATEX tested |

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-Jul-24