

1. Part Number: 8045P-1000

2. Photograph:



3. Product Information:

Description:

SteriWare PharmaScoop® PP 1000ml

Use:

Single use, disposable scoop

Material of Construction:

Polypropylene – rigid, opaque

- FDA acceptable CFR 177.1520
- EU Regulations 10/2011 compliant
- EC regulation 1935/2004 compliant
- BSE/TSE free

Method of Construction:

Single piece, injection moulded

Nominal Scoop Volume:

1000ml

Overall Length:

324mm

Handle Length:

149mm

Maximum Width:

150mm

Nominal Weight:

138g

Moulding Environment:

Class 100,000 Cleanroom

Primary Packaging Environment:

Class 10,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE/PET bag)

Method of Sterilisation:

N/A

Number of Scoops per Box:

25

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-Nov-22