

	<h1>Product Specification</h1>	Specification Number	8045P-2500
		Revision	04

Part No. 8045P-2500

Description SteriWare PharmaScoop® – 2500ml

Material of Construction: Polypropylene – rigid, opaque

- FDA acceptable CFR 177.1520
- EU Regulations 10/2011 compliant
- EC regulation 1935/2004 compliant



Method of Construction: Single piece, injection moulded

Nominal Volume: 2500ml

Overall Length: 357mm

Handle Length: 147mm

Maximum Width: 166mm

Nominal Weight: 184g

Moulding & Packing Environment: Class 100,000 Cleanroom

BSE/TSE Free Yes


Individually Bagged? Yes (heat sealed PE bag)

Method of Sterilisation: N/A

Number of Scoops per Box: 10 (Supplied as 6 boxes of 10 pcs)

Recommended Storage Conditions: Dry and ambient temperature

SteriWare and PharmaScoop are registered Trademarks belonging to Sampling Systems Ltd.

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 27-Mar-15
Page 1 of 1	Sampling Systems Ltd. 4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com	