

	<h1>Product Specification</h1>	Specification Number	8045P-500
		Revision	02

Part No. 8045P-500

Description SteriWare PharmaScoop® – 500ml

Material of Construction: Polypropylene – rigid, opaque

- FDA acceptable CFR 177.1520
- EU Regulations 10/2011 compliant
- EC regulation 1935/2004 compliant



Method of Construction: Single piece, injection moulded

Nominal Scoop Volume: 500ml

Overall Length: 249mm

Handle Length: 118mm

Maximum Width: 97mm

Nominal Weight: 71g

Moulding & Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE bag)


Method of Sterilisation: N/A

Number of Samplers per Box: 25 (Supplied as 4 boxes of 25 pcs)

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

SteriWare and PharmaScoop are registered Trademarks belonging to Sampling Systems Ltd.

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 27-Mar-15
----------------------------	--	--