

**1. Part Number:** 8049H-2000

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare® Jug 2000ml
Use:	Single use, disposable jug
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA 21 CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Internal Diameter:	129mm
Max. Diameter of Base:	157mm
Overall Height:	167mm
Nominal Weight:	232g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packaging Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
BSE/TSE Free	Yes
Number of Samplers per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature

- Shelf Life: 5 years from the date of manufacture
- ATEX: This sampler has not been ATEX tested – Contact Sampling Systems for further information
- 4. Change Control:** Sampling Systems operates a policy of continual improvement. To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)
- 5. Country of Origin:** Made and packed in the UK.
- 6. Repeated Use:** This is a single-use, disposable product. Repeated use will lead to failure.
- 7. Approval:** Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 10-Jul-24