

**1. Part Number:** 8050H-500

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare® LiquiThief - 500mm
Use:	Single use, disposable liquid sampler
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA 21 CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• Conforms to REACH regulation 1907/2006</li> <li>• BSE/TSE free</li> </ul>
	O-ring Seal – virgin Silicone Rubber <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.2600</li> </ul>
Method of Construction:	Ultrasonic welding
Overall Length	
Fully closed:	524mm
Fully extended:	885mm
Diameter of Outer Tube:	21mm
Diameter of hole in tip:	1.5mm
Width of Handle:	60mm
Nominal Weight of Sampler:	82g
Moulding Environment:	Class 100,000 Cleanroom
Assembly & Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A

Quantity per Box: 20  
Recommended Storage Conditions: Dry and ambient temperature  
Shelf Life: 5 years from the date of manufacture

#### 4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

#### 5. Country of Origin:

Made and packed in the UK.

#### 6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

#### 7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 2-Apr-25