


1. Part Number: 8051H-1100

2. Photograph:



3. Product Information:

Description:	SteriWare® Jumbo Pipette 1100mm
Use:	Single use, disposable liquid sampler
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • Conforms to FDA 21 CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free <p>A blue dye is added to one of the tips. The blue dye is certified:</p> <ul style="list-style-type: none"> • FDA acceptable for food contact • Free of materials of animal origin
Method of Construction:	Ultrasonic welding
Nominal Length	1090mm
Diameter of Outer Tube:	21mm
Diameter of hole in tip:	Handle End (Blue) = 1.5mm Product End (Natural) = 3mm
Nominal Weight of Sampler:	123g
Moulding Environment:	Class 100,000 Cleanroom
Assembly & Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Number of Samplers per Box:	37
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 11-Jul-24