


1. Part Number: 8055H-1000S

2. Photograph:



3. Product Information:

Description:	SteriWare® ViscoThief - 1000mm - Sterile
Use:	Single use, disposable liquid sampler (for high viscosity liquids)
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • Conforms to FDA 21 CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • Conforms to REACH regulation 1907/2006 • BSE/TSE free <p>O-ring Seal – virgin Silicone Rubber</p> <ul style="list-style-type: none"> • Conforms to FDA CFR 177.2600
Method of Construction:	Ultrasonic welding
Nominal Length	
Fully closed:	991mm
Fully extended:	1818mm
Diameter of Outer Tube:	21mm
Diameter of hole in tip:	10mm
Width of Handle:	103mm
Nominal Weight of Sampler:	158g
Moulding Environment:	Class 100,000 Cleanroom
Assembly & Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25kGy minimum)
Number of Samplers per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 11-Jul-24