



**1. Part Number:** 8055H-1000

**2. Photograph:**



**3. Product Information:**

<b>Description:</b>	SteriWare® ViscoThief - 1000mm
<b>Use:</b>	Single use, disposable liquid sampler (for high viscosity liquids)
<b>Material of Construction:</b>	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA 21 CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• Conforms to REACH regulation 1907/2006</li> <li>• BSE/TSE free</li> </ul> <p>O-ring Seal – virgin Silicone Rubber</p> <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.2600</li> </ul>
<b>Method of Construction:</b>	Ultrasonic welding
<b>Nominal Length</b>	
Fully closed:	991mm
Fully extended:	1818mm
<b>Diameter of Outer Tube:</b>	21mm
<b>Diameter of hole in tip:</b>	10mm
<b>Width of Handle:</b>	103mm
<b>Nominal Weight of Sampler:</b>	158g
<b>Moulding Environment:</b>	Class 100,000 Cleanroom
<b>Assembly &amp; Primary Packing Environment:</b>	Class 10,000 Cleanroom
<b>Individually Bagged?</b>	Yes (heat sealed PE/PET bag)
<b>Method of Sterilisation:</b>	N/A
<b>Number of Samplers per Box:</b>	20
<b>Recommended Storage Conditions:</b>	Dry and ambient temperature
<b>Shelf Life:</b>	5 years from the date of manufacture
<b>ATEX:</b>	This product has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 11-Jul-24