


**1. Part Number:** 8060H-1520S

**2. Photograph:**



**3. Product Information:**

**Description:** SteriWare PowderThief - 1500mm – 20ml Tip - Sterile

**Use:** Single use, disposable powder sampler

**Material of Construction:** HDPE (High Density Polyethylene) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

**Method of Construction:** Ultrasonic welding

**Nominal Length:** 1500mm

**Diameter of Outer Tube:** 21mm

**Nominal Tip Size:** 20ml

**Width of Handle:** 103mm

**Nominal Weight of Sampler:** 238g

**Moulding Environment:** Class 100,000 Cleanroom

**Assembly & Primary Packing Environment:** Class 10,000 Cleanroom

**Individually Bagged?** Yes (heat sealed PE/PET bag)

**Method of Sterilisation:** Gamma Irradiation (25-45kGy)

**Number of Samplers per Box:** 20

**BSE/TSE Free:** Yes

**Recommended Storage Conditions:** Dry and ambient temperature

**Shelf Life:** 5 years from the date of manufacture

**ATEX:** This device has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 11-Jul-24