

1. Part Number: 8060H-501S

2. Photograph:



3. Product Information:

Description:	SteriWare PowderThief - 500mm – 1ml Tip – Sterile
Use:	Single use, disposable powder sampler
Material of Construction:	HDPE (High Density Polyethylene) - virgin <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Ultrasonic welding
Nominal Length:	511mm
Diameter of Outer Tube:	21mm
Nominal Tip Size:	1ml
Width of Handle:	103mm
Length 'L':	5mm
Nominal Weight of Sampler:	93g
Moulding Environment:	Class 100,000 Cleanroom
Assembly & Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (minimum dose 25 kGy)
Number of Samplers per Box:	20
BSE/TSE Free:	Yes
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 11-Jul-24