

1. Part Number: 8060H-505

2. Photograph:



3. Product Information:

| | |
|---|--|
| Description: | SteriWare PowderThief - 500mm – 5ml Tip |
| Use: | Single use, disposable powder sampler |
| Material of Construction: | HDPE (High Density Polyethylene) - virgin <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free |
| Method of Construction: | Ultrasonic welding |
| Nominal Length: | 511mm |
| Diameter of Outer Tube: | 21mm |
| Nominal Tip Size: | 5ml |
| Width of Handle: | 103mm |
| Length 'L': | 31mm |
| Nominal Weight of Sampler: | 93g |
| Moulding Environment: | Class 100,000 Cleanroom |
| Assembly & Primary Packing Environment: | Class 10,000 Cleanroom |
| Individually Bagged? | Yes (heat sealed PE/PET bag) |
| Method of Sterilisation: | N/A |
| Number of Samplers per Box: | 20 |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Shelf Life: | 5 years from the date of manufacture |
| ATEX: | This device has not been ATEX tested |

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-Nov-22