

1. Part Number: 8064H-1000S

2. Photograph:



3. Product Information:

Description:	SteriWare® Powder Probe 1000mm - Sterile
Use:	Single use powder sampler
Material of Construction:	High Density Polyethylene (HDPE) – virgin <ul style="list-style-type: none"> • Conforms to FDA 21 CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Ultrasonic welding
Nominal Length:	945mm
Blade Length:	845mm
Diameter of Blade:	22mm
Width of opening in blade:	10mm
Nominal Weight of Sampler:	123g
Moulding Environment:	Class 100,000 Cleanroom
Welding & Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Quantity per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Jul-24