


1. Part Number: 8066A-1

2. Photograph:




Embossed information:

ROWS	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
TABLETS	1	3	6	10	15	21	28	36	45	55	66	78	91	105	120

3. Product Information:

Description:	SteriWare® Tablet Counting Triangle
Use:	Single use tablet tray for counting round tablets
Material of Construction:	High Impact Polystyrene (HIPS) – rigid, white <ul style="list-style-type: none"> • Conforms to FDA 21 CFR 177.1640 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Length:	230mm
Height:	13mm
Nominal Weight:	57g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Quantity per Box:	30
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested



Sampling Systems Ltd.

2 Trillennium Point, Gorse Lane, Coleshill, B46 1JU. UK.
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Jul-24