

**1. Part Number:** 8069H-500S

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare® Mono Slot Sampler – 500mm – Sterile
Use:	Single use, disposable powder sampler
Material of Construction:	HDPE (High Density Polyethylene) - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Ultrasonic welding
Nominal Length:	620mm
Diameter of Blade:	22mm
Diameter of Outer Tube:	26mm
Number of Sampling Slots:	1
Slot Positions:	



Sampling Slot Width:	8mm
Sampling Slot Length:	60mm
Nominal Weight of Sampler:	167g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25 – 50 kGy)
Quantity per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Sep-20