


1. Part Number: 8073H-300S

2. Photograph:



3. Product Information:

Description:

SteriWare Slot Sampler - 300mm - Sterile

Use:

Single use, disposable powder sampler

Material of Construction:

High Density Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction:

Ultrasonic welding

Overall Length:

420mm

Diameter of Outer Tube:

26mm


Diameter of Inner Blade:

22mm

Nominal Weight of Sampler:

100g

Slot Positions:



Number of Sampling Slots:

3

Sampling Slot Width:

11mm

Assembly & Primary Packing Environment: Class 100,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE/PET bag)

Method of Sterilisation:

Gamma Irradiation (Minimum Dose 25kGy)

Number of Samplers per Box:

20

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Jul-24