


1. Part Number: 8073H-500S

2. Photograph:



3. Product Information:

Description:	SteriWare Slot Sampler - 500mm - Sterile
Use:	Single use, disposable powder sampler
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • Conforms to FDA 21 CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Ultrasonic welding
Overall Length:	635mm
Diameter of Outer Tube:	26mm
Diameter of Inner Blade:	22mm
Nominal Weight of Sampler:	160g
Slot Positions:	



Number of Sampling Slots:	3
Sampling Slot Width:	11mm
Manufacturing Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (Minimum Dose 25kGy)
Number of Samplers per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 11-Oct-24