

**1. Part Number:** 8073J-1000S

**2. Photograph:**



**3. Product Information:**

**Description:** SteriWare Jumbo Slot Sampler - 1000mm - Sterile

**Use:** Single use, disposable powder sampler

**Material of Construction:** High Density Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

**Method of Construction:** Ultrasonic welding

**Overall Length:** 986mm

**Diameter of Outer Tube:** 40mm

**Diameter of Inner Tube:** 35mm


**Nominal Weight of Sampler:** 297g

**Number of Sampling Slots:** 3

**Sampling Slot Width:** 8mm

**Sampling Slot Length:** 60mm

**Sampling Slot Positions:**



**Assembly & Primary Packing Environment:** Class 100,000 Cleanroom

**Individually Bagged?** Yes (heat sealed PE/PET bag)

**Method of Sterilisation:** Gamma Irradiation (25-50 kGy)



**Sampling Systems Ltd.**

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.

Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

Number of Samplers per Box: 20  
Recommended Storage Conditions: Dry and ambient temperature  
Shelf Life: 5 years from the date of manufacture

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Jun-21