

1. Part Number: 8073J-1000S

2. Photograph:



3. Product Information:

Description: SteriWare Jumbo Slot Sampler - 1000mm - Sterile

Use: Single use, disposable powder sampler

Material of Construction: High Density Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction: Ultrasonic welding

Overall Length: 986mm

Diameter of Outer Tube: 40mm

Diameter of Inner Tube: 35mm


Nominal Weight of Sampler: 297g

Number of Sampling Slots: 3

Sampling Slot Width: 8mm

Sampling Slot Length: 60mm

Sampling Slot Positions:



Assembly & Primary Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (25-50 kGy)

Number of Samplers per Box: 20
Recommended Storage Conditions: Dry and ambient temperature
Shelf Life: 5 years from the date of manufacture

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Jul-24