

Specification Sheet

Specification Number: 8079A-450
Revision: 02
Date Issued: 14-Mar-24

1. Part Number: 8079A-450

2. Product Information:

Description: Sterile Bag Stand-Up 450ml

Use: Single use, disposable sampling bags

Material of Construction: Thermoplastic laminate

• Conforms to FDA 21 CFR 177.1520

• Conforms to EU 10/2011

• BSE/TSE free

Closure Type: Attachment with 2 round wires

Bag Dimensions: 114 x 229mm

Thickness of Bag 89 microns

Nominal Bag Volume: 400ml

Method of Sterilisation: Sterility is provided by irradiation. The sterilization dose ensuring a

sterility assurance level (SAL) of 10⁻³ has been established

according to ISO 11137.

Pyrogens: This product is non-pyrogenic at the endotoxin limit of 2.15

EU/device. Nonpyrogenicity is supported by endotoxin testing of randomly selected samples using the Limulus amebocyte lysate (LAL) gel assay according USP-NF <85> and <161> guidelines.

Packaging: 1000 bags per case (4 x 250)
Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years

ATEX: This product has not been ATEX tested

3. Country of Origin: Made and packed in Canada.

4. Repeated Use: This is a single-use, disposable product and should not be re-

used.

5. Approval: Signed: Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Mar-24

Note: The information given in this document has been given to us by our suppliers and is passed on in good faith, we believe the information to be correct.