

1. Part Number: 8079A-600

2. Product Information:

Description: Sterile Bag Stand-Up 600ml

Use: Single use, disposable sampling bags

Material of Construction: Thermoplastic laminate

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU 10/2011
- BSE/TSE free

Closure Type: Attachment with 2 round wires

Bag Dimensions: 152 x 229mm

Thickness of Bag: 89 microns

Nominal Bag Volume: 600ml

Method of Sterilisation: Sterility is provided by irradiation. The sterilization dose ensuring a sterility assurance level (SAL) of 10^{-3} has been established according to ISO 11137.

Pyrogens: This product is non-pyrogenic at the endotoxin limit of 2.15 EU/device. Nonpyrogenicity is supported by endotoxin testing of randomly selected samples using the Limulus amoebocyte lysate (LAL) gel assay according to USP-NF <85> and <161> guidelines.

Packaging: 1000 bags per case (4 x 250)

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years

ATEX: This product has not been ATEX tested

3. Country of Origin: Made and packed in Canada.

4. Repeated Use: This is a single-use, disposable product and should not be re-used.

5. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Mar-24

Note: The information given in this document has been given to us by our suppliers and is passed on in good faith, we believe the information to be correct.