

1. Part Number: 8081P-3000S

2. Photograph:



3. Product Information:

Description:	SteriWare Tray 3000ml - sterile
Use:	Single use, disposable tray
Material of Construction:	Polypropylene (PP) - rigid <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulation 10/2011 • Conforms to EC Regulation 1935/2004 • BSE/TSE free
Method of Construction:	Single Piece, injection moulded
Overall Dimensions:	318 x 267 x 55mm
Nominal Volume:	3000ml
Nominal Weight:	180g
Assembly & Primary Packing Environment:	Moulded and packed in a class 100,000 cleanroom.
Individually Bagged?	Yes (heat sealed PE bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Number of Trays per Box:	30
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	2 years from the date of manufacture
ATEX:	This sampler has not been ATEX tested – Contact Sampling Systems for further information

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable product. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Oct-20