

**1. Part Number:** 8081P-3000

**2. Photograph:**



### 3. Product Information:

<b>Description:</b>	SteriWare Tray 3000ml
<b>Use:</b>	Single use, disposable tray
<b>Material of Construction:</b>	Polypropylene (PP) - rigid <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulation 10/2011</li> <li>• Conforms to EC Regulation 1935/2004</li> <li>• BSE/TSE free</li> </ul>
<b>Method of Construction:</b>	Single Piece, injection moulded
<b>Overall Dimensions:</b>	318 x 267 x 55mm
<b>Nominal Volume:</b>	3000ml
<b>Nominal Weight:</b>	180g
<b>Assembly &amp; Primary Packing Environment:</b>	Moulded and packed in a class 100,000 cleanroom.
<b>Individually Bagged?</b>	Yes (heat sealed PE bag)
<b>Method of Sterilisation:</b>	N/A
<b>Number of Trays per Box:</b>	40
<b>Recommended Storage Conditions:</b>	Dry and ambient temperature
<b>Shelf Life:</b>	2 years from the date of manufacture
<b>ATEX:</b>	This sampler has not been ATEX tested – Contact Sampling Systems for further information

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable product. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Oct-20