


**1. Part Number:** 8088A-03

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare® Sample Spoon 3ml
Use:	Single use, disposable spoon
Material of Construction:	High Impact Polystyrene (HIPS) – rigid, white <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1640</li> <li>• Conforms to EU regulation 10/2011</li> <li>• Conforms to EC regulation 1935/2004</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Nominal Bowl Volume:	3ml
Length:	173mm
Bowl Length:	36mm
Bowl Width:	25mm
Nominal Weight:	4g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Number of Samplers per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

[change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable device. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Jul-24