

**1. Part Number:** 8088B-10S

**2. Photograph:**



### 3. Product Information:

|                                 |  |
|---------------------------------|--|
| Description:                    | SteriWare Powder Spoon Blue - 10ml - Sterile   |
| Use:                            | Single use, disposable spoon   |
| Material of Construction:       | Polypropylene (PP) – rigid, blue <ul style="list-style-type: none"> <li>• FDA acceptable CFR 177.1520</li> <li>• EU Regulations 10/2011 compliant</li> <li>• EC Regulations 1935/2004 compliant</li> <li>• BSE/TSE free</li> </ul> |
| Method of Construction:         | Single piece, injection moulded  |
| Nominal Bowl Volume:            | 10ml   |
| Length:                         | 185mm  |
| Bowl Length:                    | 48mm   |
| Bowl Width:                     | 35mm   |
| Nominal Weight:                 | 5g   |
| Moulding Environment:           | Class 100,000 Cleanroom  |
| Primary Packing Environment:    | Class 10,000 Cleanroom   |
| Individually Bagged?            | Yes (heat sealed PE/PET bag)   |
| Method of Sterilisation:        | Gamma Irradiation (25-45kGy)   |
| Number of Samplers per Box:     | 100  |
| Recommended Storage Conditions: | Dry and ambient temperature  |
| Shelf Life:                     | 5 years from the date of manufacture   |
| ATEX:                           | This product has not been ATEX tested  |

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Jul-24