


	<h1>Product Specification</h1>	Specification Number	8110A-120
		Revision	01

Part No. 8110A-120
Description Aluminium Bottle with tamper evident cap
Nominal Volume: 120ml



Material of Construction: Bottle: Aluminium
Cap: Polypropylene with aluminium seal

Bottle Height with Cap: 97mm
Bottle Height: 91mm
Bottle Diameter: 57mm
Internal Diameter of Neck: 25mm
Nominal Weight: 45g (Bottle & Cap)
Recommended Storage Conditions: Dry and ambient temperature

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 11-Apr-13
Page 1 of 1	Sampling Systems Ltd. 4-5 Forge Mills Park, Station Road, Coleshill, Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com	