

**1. Part Number:** 8111H-100S

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare Powder Funnel 100mm - Sterile
Use:	Single use, disposable powder spatula & cover
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> <li>• FDA acceptable CFR 177.1520</li> <li>• EU Regulations 10/2011 compliant</li> <li>• EC Regulations 1935/2004 compliant</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Max. Diameter:	117mm
ID at Top of Funnel:	100mm
Overall Height:	88mm
OD at Bottom of Spout:	23mm
ID at Bottom of Spout:	20mm
Nominal Weight:	19g
Moulding Environment:	Class 100,000 Cleanroom
Packing Environment	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Quantity per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture

ATEX:

This product has not been ATEX rated – Contact Sampling Systems for details

#### 4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

#### 5. Country of Origin:

Made and packed in the UK.

#### 6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

#### 7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-Nov-22