

1. Part Number: 8111H-200

2. Photograph:



3. Product Information:

Description:	SteriWare Powder Funnel 200mm
Use:	Single use, disposable powder spatula & cover
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • FDA acceptable CFR 177.1520 • EU Regulations 10/2011 compliant • EC Regulations 1935/2004 compliant • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Max. Diameter:	200mm
Overall Height:	165mm
OD at Bottom of Spout:	43mm
ID at Bottom of Spout:	40mm
Nominal Weight:	86g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A

Quantity per Box: 10
Recommended Storage Conditions: Dry and ambient temperature
Shelf Life: 5 years from the date of manufacture
ATEX: This product has not been ATEX rated

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 23-Jul-24